

If you would like to significantly reduce the amount of time you will need to wear braces or aligners and reduce the pain usually experienced when adjustments to your braces are made or new aligners are started, then OrthoAlign is what you need to use along with your orthodontic treatment.


Helps move teeth up to 50% faster for anyone wearing braces! It is a comfortable, safe and lightweight device.

50%



Accelerate your smile journey with Orthodontic Accelerator faster, more comfortable, and effective treatment for a confident smile!

Another benefit of accelerated orthodontics is that any type of braces can be used. This includes metal braces (the most basic and widely used type), ceramic braces (white or clear brackets that are less visible), and lingual braces (placed behind the teeth so they cannot be seen).



Perfect Smile

Orthodontics now offers this incredible technology to help ease any discomfort and speed your orthodontic treatment.

CONTACT CENTURY MEDICAL, INC.

-  1712 Pioneer Ave. Ste 500
Cheyenne, WY 82001
USA
-  +1(307)-222-9655
-  info@cmedicali.com
-  www.cmedicali.com

Pulse vibration device

TRIFOLD BROCHURE TEMPLATE

*Get your beautiful,
healthy smile – faster!*

Device emits tiny, gentle vibrations, known as micropulses.
The micropulses accelerate the changes in your bone that allow your teeth to move
(also known as "bone remodeling").



Accelerate
tooth movement




Reduce pain during
orthodontic treatment



Align teeth and
braces accurately



Intensive fixed
orthodontic effect



What allows orthodontics to move teeth is the careful manipulation of force that is used to guide the teeth into a new, improved position and better equilibrium.

CASE

Midline Correction
6.5 months



Anterior open bite
7 months



Anterior Crossbite
10 months



Upper & Lower Crowding
5.5 months



You wear it just 5~20 minutes a day, every day to speed up tooth movement.

Helps you achieve more
dramatic results in less time